

Agilent 1100 Binary Pump



Agilent 1100 HPLC Binary pump comes with new pistons, seals and fully refurbished and tested.

Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Manufacturer: [Agilent](#)

Description

Agilent 1100 Binary Pump

The binary pump comprises two identical pumps integrated into one housing. It provides gradient generation by high-pressure mixing.

Specifications:

- Flow Range: .001-5 ml/min
- Flow Precision: less than 0.3 % RSD
- Pressure Range: 0-400 bar (0-5880 psi)
- Pressure Pulsation Less than 2% amplitude
- Recommended PH Range: 1.0-12.5

Includes:

- CAN Cable
- Manual
- **NINETY (90) DAY WARRANTY**

One (1) Year Warranty Available

