

Agilent 1100 Column Compartment



This Agilent / HP 1100 Series Column Compartment comes fully refurbished and tested.

Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Manufacturer: [Agilent](#)

Description

Agilent / HP 1100 Series Column Compartment

Features:

- Peltier heating and cooling from 10 degrees below ambient up to 80 degrees Centigrade.
- Holds up to three 30-cm columns
- Two independently programmable heat exchangers.
- Electronic column-identification module.

Specifications:

- Temperature range: 10 degrees below ambient to 80 degrees Centigrade
- Temperature stability + or – 0.15 degrees centigrade
- Temperature accuracy: + or – 0.8 degrees centigrade
- Column Capacity: Three 30 cm columns

Includes:

- CAN Cable
- Manual

- **NINETY (90) DAY WARRANTY**

One (1) Year Warranty Available