

Agilent 1100 Diode Array Detector



Agilent HPLC 1100 Diode Array Detector. It comes fully refurbished and tested.

Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Manufacturer: [Agilent/HP](#)

Description

Agilent 1100 Diode Array Detector

Features:

- Expanded Detection Range through deuterium and tungsten lamps.
- 1024 diodes and a 1 nm slit for high spectral resolution
- A programmable slit range of 1 to 16 nm

Specifications:

- Detector type: 1024 element diode array
- Short term noise: + and – 1.0×10^{-5} AU at 254 and 750 nm
- Drift: 2×10^{-3} AU/h at 254 nm
- Linearity: Greater than 2 AU upper limit
- Wavelength range 190 to 950 nm
- Wavelength Accuracy: + and – 1 nm
- Flow Cells: Standard

Includes:

- CAN Cable
- Manual
- **NINETY (90) DAY WARRANTY**

One (1) Year Warranty Available