

Agilent 1100 Quaternary Pump



This Agilent / HP 1100 Series Quaternary Pump comes fully refurbished and tested.

Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Manufacturer: [Agilent](#)

Description

Agilent / HP 1100 Series Quaternary Pump

Specifications:

- Hydraulic system – Dual Plunger in series pump
- Flow Range – 0.2-10.0 ml/min
- Flow Precision – <0.3% RSD at 1 ml/min
- Pressure 0-40 MP (0-5880 psi) up to 5 ml/min
- Pressure 0-20 MP (0-2950 psi) up to 10 ml/min
- Pressure pulsation < 2% amplitude at 1 ml/min isopropanol
- **NINETY (90) DAY WARRANTY**

Includes:

- CAN Cable
- Manual

ONE (1) YEAR WARRANTY Available

