

Agilent 1100 Series Autosampler G1313



This Agilent 1100 Series Autosampler with 100 position Sample Tray comes fully refurbished and tested with 1-Year Warranty.

Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Sales price with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Manufacturer: [Agilent](#)

Description

Agilent 1100 Series Autosampler G1313

Features:

- Communications: Controller Area Network (CAN), GPIB (IEEE-448), Serial (RS232C)
- GLP features: Early maintenance feedback (EMF), records of maintenance and errors
- Safety features: Leak detection, safe leak handling, and error detection and display
- Injection range: From 0.1 – 100 µl in 0.1 µl increments
- Replicate injections: From 1 – 99 from one vial
- Injection cycle time: Typically 50 s, depending on draw speed and volume
- Sample capacity: 100 x 2 ml vials in tray

Includes:

- Sample tray
- CAN Cable
- Manual
- **NINETY (90) DAY WARRANTY**

ONE (1) YEAR WARRANTY Available

