

## Agilent 1100 Series Isocratic Pump G1310



This Agilent / HP 1100 Series Isocratic Pump (G1310) comes fully refurbished and tested.

Rating: Not Rated Yet

**Price:**

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Manufacturer: [Agilent](#)

Description

### Agilent 1100 Series Isocratic Pump G1310

#### Specifications:

- Hydraulic system: Dual Piston in series pump
- Flow Range: 0.05 – 10.0 ml/min
- Flow Precision: <=0.07% RSD
- Pressure: 0-40 MP (0-5880 psi) up to 5 ml/min
- Pressure: 0-20 MP (0-2950 psi) up to 10 ml/min
- Pressure pulsation: < 2% amplitude at 1 ml/min isopropanol

#### Includes:

- Isocratic Pump
- CAN Cable
- Manual
- **NINETY (90) DAY WARRANTY**

**One (1) Year Warranty Available**

