

Agilent 1100 Variable Wavelength Detector



HP/Agilent 1100 Series Variable Wavelength Detector (VWD) model G1314A comes fully refurbished and tested.

Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Manufacturer: [Agilent](#)

Description

HP/Agilent 1100 Variable Wavelength Detector VWD

Features:

- Wavelength Range: 190-600 nm
- Wavelength Accuracy: +/- 1 nm
- Band width: 6.5 nm
- Double beam photometer with Deuterium lamp
- Linearity: >2 AU (5%) upper limit
- Drift: 3×10^{-4} AU/hr at 254 nm
- Analog output connection

Includes:

- Standard flow cell
- CAN cable
- **NINETY (90) DAY WARRANTY**

ONE (1) YEAR WARRANTY Available

