

Agilent 1200 Series ELSD Detector



This Agilent 1200 Series ELSD Detector comes fully refurbished and tested.

Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Manufacturer: [Agilent](#)

Description

Agilent 1200 Series ELSD Detector

Features:

- Detection: High Sensitivity Photomultiplier
- Light Source: Selected high Efficiency Blue LED
- Temperature Range: Ambient to 100 degrees C
- Gas Flow Control: Computer controlled nebulization gas flow and auxiliary gas flow
- Standard Flow Nebulizer: 500 ul/min to 2 ml/min
- **NINETY (90) DAY WARRANTY**
- Satisfaction Guaranteed!

One (1) Year Warranty Available

