

## Agilent 1200 Series High Performance Autosampler SL G1367D



This Agilent 1200 Series High Performance Autosampler SL G1367D comes fully tested.

Rating: Not Rated Yet

**Price:**

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Manufacturer: [Agilent](#)

Description

### Agilent 1200 Series High Performance Autosampler SL G1367D

**Features:**

- **NEW – Has not been used!**
- Early maintenance feedback (EMF)
- Controller-area network (CAN). RS232C, APG-remote standard
- Leak detection and safe leak handling
- Injection Range: 0.1 – 40 µl in 0.1 µl increments Up to 100 µl with extended injection
- Precision: Typically < 0.25% RSD from 5 – 40 µl, Typically < .5% RSD
- Injection cycle time: Typically < 17 s using standard conditions
- Includes well plate tray, CAN cable and Manual
- **NINETY (90) DAY WARRANTY**
- Satisfaction Guaranteed!

