

Agilent 1260 Isocratic Pump G1310B



This Agilent 1260 Isocratic Pump G1310B comes fully tested.

Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Manufacturer: [Agilent](#)

Description

Agilent 1260 Isocratic Pump G1310B

Features:

- **NEW – Has not been used!**
- Dual piston in series pump with servo-controlled variable stroke drive, floating pistons
- Set points 0.001 – 5 mL/min, in 0.001 mL/min increments
- Flow range 0.2 – 10.0 mL/min
- Flow precision: less than 0.07 % RSD or less than 0.02 min SD, whatever is greater
- Flow accuracy $\pm 1\%$ or 10 $\mu\text{L}/\text{min}$ whatever is greater
- Pressure operating range: Operating range up to 60 MPa (600 bar, 8700 psi) up to 5 mL/min Operating range up to 20 MPa (200 bar, 2950 psi) up to 10 mL/min
- Includes CAN Cable and Manual
- **NINETY (90) DAY WARRANTY**

