

HP 1050 Quaternary HPLC Pump



The HP 1050 Quaternary HPLC Pump is in excellent condition.

Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Manufacturer: [Agilent/HP](#)

Description

HP 1050 Quaternary HPLC Pump

Features:

- Refurbished and Tested at GenTech.
- HPIB Communications available.
- Includes **NINETY (90) DAY WARRANTY.**

One (1) Year Warranty Available

Specifications:

- Flow range from .001 to 9.999 ml per minute
- Pressure range from 0 to 400 bar up to 5 ml per minute
- PH range of 2.3 to 12.5
- Composition range from 0 to 100% in .1% steps
- Battery backed storage of up to 10 methods

