

LabAlliance Prep 24 or Prep 100 Preparative Pump



The Prep Pumps have been specifically designed for Preparative Chromatography applications. The units provide high flow rates and high pressure capabilities, with consistent performance to minimize the potential for column damage.

Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Manufacturer: [LabAlliance](#)

Description

The Prep Pumps have been specifically designed for **Preparative Chromatography applications**. The units provide high flow rates and high pressure capabilities, with consistent performance to minimize the potential for column damage.

The units deliver a maximum flow rate of 24 or 100 mL/min. at pressures to **6000 psi**.

Standard controls include RS-232 serial communication port, Run/Stop inputs and remote frequency & voltage controls. Real-time pressure monitoring with user selectable upper & lower pressure limits is standard. The Prep pumps offer easy solvent changeover and priming capability. Piston wash and simple seal change keep maintenance costs down.

Includes **THREE (3) YEAR WARRANTY**.

Available in PEEK or Stainless Steel Tubing.

P40PFT01 PEEK Prep 100 Preparative Pump

P40SFT01 Stainless Prep 100 Preparative Pump

P10PFT01 PEEK Prep 24 Preparative Pump

P10SFT01 Stainless Prep 24 Preparative Pump

