

LabAlliance Series II Digital Pump



The Series II pump consists of a series of modular sub-assemblies packaged to meet basic metering and chromatography application requirements. The extensive use of digital logic and a micro-stepped motor allow for precise and flexible control. The user interface consists of an LED display and keypad. The use of the self-flush head and a three-year warranty ensure years of trouble-free service. Method development is easy with Autoprime™ purging. The Series II gives the modern chromatographer solid, reliable performance at a value price. When combined with a Series III pump, keypad controller, detector, and injector, you have an economical, flexible binary gradient system.

Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Manufacturer: [LabAlliance](#)

Description

The **Series II pump** consists of a series of modular sub-assemblies packaged to meet basic metering and chromatography application requirements. The extensive use of digital logic and a micro-stepped motor allow for **precise and flexible control**. The user interface consists of an LED display and keypad. The use of the self-flush head and a three-year warranty ensure years of trouble-free service. Method development is easy with Autoprime™ purging. The Series II gives the modern chromatographer solid, reliable performance at a value price. When combined with a Series III pump, keypad controller, detector, and injector, you have an economical, flexible binary gradient system.

Includes **THREE (3) YEAR WARRANTY**.

Available in Bioclean or Stainless Steel Tubing.

210PFP01 PEEK Series II Digital Pump, 0.01-9.99 mL/min.

210SFP01 Stainless Series II Digital Pump, 0.01-9.99 mL/min.