

LabAlliance Series III Digital Pump



The Series III pump consists of a series of modular sub-assemblies packaged to meet any isocratic application requirements. The extensive use of digital logic and a micro-stepped motor allow for precise and flexible control. The user interface consists of an LED display and keypad. The self-flush head and a three-year warranty ensure years of trouble-free service. Method development is easy with Autoprime™ purging. Upper and lower pressure limits can be easily set, and the pressure can be monitored on the LED display. External control gives flexible performance at an affordable price.

Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Manufacturer: [LabAlliance](#)

Description

The **Series III pump** consists of a series of modular sub-assemblies packaged to **meet any isocratic application requirements**. The extensive use of digital logic and a micro-stepped motor allow for precise and flexible control. The user interface consists of an LED display and keypad. The self-flush head and a **three-year warranty ensure years of trouble-free service**. Method development is easy with Autoprime™ purging. Upper and lower pressure limits can be easily set, and the pressure can be monitored on the LED display. External control gives flexible performance at an affordable price.

Includes **THREE (3) YEAR WARRANTY**.

Available in Bioclean or Stainless Steel Tubing.

310PFT01 PEEK Series III Digital Pump, 0.01-9.99 mL/min.

310SFT01 Stainless Series III Digital Pump, 0.01-9.99 mL/min.

340PFT01 PEEK Series III Digital Pump, 0.04-40 mL/min.

340SFT01 Stainless Series III Digital Pump, 0.04-40 mL/min.

305PFT01 PEEK Series III Digital Pump, 0.001-5 mL/min.

305SFT01 Stainless Series III Digital Pump, 0.001-5 mL/min.