

New Leybold TurboVac 151 Turbo pump



WC014709 PUMP TMP 151 VACUUM TURBO PUMP

Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Manufacturer: [Leybold](#)

Description

New Leybold TurboVac 151 Turbo Pump

Does your mass spec need a new turbo pump? Do you need a spare one?

GenTech Scientific has lots of Brand New Turbo Pumps in stock! Give us a call!

The Turbovac 151 pumps for ABI Maldi TOF mass spec's:

- Ultimate pressure 10-10 mbar
- Speed: 36,000 rpm
- Forevacuum port: DN 40 KF
- Purge gas/vent port: DN 10 KF
- Cooling Water Connection: 10mm
- Weight: 12 kg
- **NINETY (90) DAY WARRANTY**