

Precision Differential Refractive Index Detector IOTA 2



High stability in RI detection The unique and patented system provides a fully symmetrical in line beam, specifically designed to avoid all external vibration problems.

Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Manufacturer: [LabAlliance](#)

Description

High stability in RI detection The unique and patented system provides a fully symmetrical in line beam, specifically designed to avoid all external vibration problems.

Fully automatic Real auto-zero system, microprocessor controlled, push-button reference cell filling, RS232.

For analytical or preparative chromatography PRECISION Instruments IOTA 2 Differential Refractometer is designed to offer high detection stability and sensitivity. Three different cells are available for analytical or preparative chromatography.

ONE (1) YEAR WARRANTY.