

## Sensivate PCR II Post-Column Reactor



This Sensivate PCR II Post-Column Reactor system which includes two series I pump and Dual Reactor is totally modular. This ensures fast and easy set up, maintenance, and repair. The “All-In-One” Integrated Fluidics unit is a combination Pulse Damper, Prime Purge Valve, Outlet Filter, and Transducer. This eliminates tubing and fittings between components and makes for fewer leak points and very low dead volume. All fittings are conveniently positioned in the front of the units. Add post-column capability to any LC system. Dual Chamber, 3 configurations offered, Heated, See below

Rating: Not Rated Yet

**Price:**

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Manufacturer: [LabAlliance](#)

### Description

This **Sensivate PCR II Post-Column Reactor** system which includes two series I pump and Dual Reactor is totally modular. This ensures fast and easy set up, maintenance, and repair. The “All-In-One” Integrated Fluidics unit is a combination Pulse Damper, Prime Purge Valve, Outlet Filter, and Transducer. This eliminates tubing and fittings between components and makes for fewer leak points and very low dead volume. All fittings are conveniently positioned in the front of the units. Add **post-column capability to any LC system**. Dual Chamber, 3 configurations offered, Heated, See below

Includes **ONE (1) YEAR WARRANTY**

**PCR2-R050** Dual Chamber, .5 ml Heated

**PCR2-R050-R015** Dual Chamber, .5 ml Heated, .15 ml Ambient

**PCR2-R140-R140** Dual Chamber, 1.4 ml Heated, 1.4 ml Ambient